


The EXP2 Pre-Column and In-Line Hand-Tight Filters are ideal for protecting small-particle UHPLC columns and systems at extreme pressures. Such techniques analyze samples in the most demanding applications, decreasing the life of expensive hardware. EXP2 Filters help extend the life of a UHPLC system and protect your investment without sacrificing performance. Free-Turn<sup>®</sup> architecture allows the user to change cartridges by hand without breaking fluid connections on the holder inlet/outlet.

EXP2 Pre-Column Filter holders are available with Titanium Hybrid Ferrules for easy direct connection to any 10-32 port. The filter holder is a complete package and includes fittings to provide repeated zero-dead-volume connections.

## Advantages of EXP<sup>®</sup>2 Filters

- Approved for use at 20,000+ psi (1,400+ bar)
- Hand-tight filter replacement - no tools needed
- Auto-adjusting ZDV connections
- Depth filtration maximizes capacity
- Low-volume, low-dispersion cartridges
- Hardened stainless steel end cap eliminates galling
- New design offers more cost efficient cartridges



## EXP<sup>®</sup>2 Filter Holders

- 15-04-05242** EXP2 Pre-Column Filter Holder
- 15-04-05258** EXP2 In-Line Filter Holder

The above holders include fittings

## EXP<sup>®</sup>2 Filter Cartridges

5 PACK	10 PACK	
<b>15-04-05263</b>	<b>15-04-05264</b>	0.2µm, 3mm Cartridge
<b>15-04-05261</b>	<b>15-04-05262</b>	0.5µm, 3mm Cartridge
<b>15-04-05259</b>	<b>15-04-05260</b>	2µm, 3mm Cartridge

## EXP<sup>®</sup> Filter Cartridges (Old Style)

5 PACK	10 PACK	
<b>15-04-03096</b>	<b>15-04-03097</b>	0.2µm Cartridge
<b>15-04-03093</b>	<b>15-04-03094</b>	0.5µm Cartridge

